

Directions for use

ZEDI*X*CITE 330/418 FLUOROGENIC TG2-ASSAY KIT

Isopeptidase-activity based TG2 assay
(Fluorogenic: $\lambda_{\text{ex}} = 330 \text{ nm}$; $\lambda_{\text{em}} = 418 \text{ nm}$)

Art.-No. F014

For Research & Development Only

Revision Number: 2.1

Release Date: 2023-12-08



Zedira GmbH, Roesslerstrasse 83, 64293 Darmstadt, Germany
www.zedira.com contact@zedira.com



Watch the video on YouTube, describing
how to run ZEDXCITE FLUOROGENIC TG2-ASSAY KIT F014

Table of contents

1. Assay principle	3
2. Intended use	3
3. Test sample	3
4. Reagents in the kit.....	5
5. Reagent preparation, storage, and stability	5
6. Schematic assay overview	7
7. Procedure and Equipment	9
8. Number of sample measurements.....	10
9. Results	10
10. Calibration.....	12
11. Reference Range.....	14
12. Kinetic parameters	14
13. Limitations.....	14
14. Precision	15
References	15

1. Assay principle

TISSUE TRANSGLUTAMINASE (TG2, transglutaminase 2) is a multifunctional enzyme best known for its crosslinking activity, but it can also cleave isopeptide bonds. This feature is used to provide an **easy to handle, robust, sensitive, and precise fluorogenic assay** to measure TG2 activity. The assay is suitable for research purposes, drug discovery programs, and quality assurance.

TISSUE TRANSGLUTAMINASE cleaves a dark quencher molecule from the side chain of a TG2 specific peptide (SUBSTRATE REAGENT **SR**). The TG2-substrate complex is cleaved by the incorporation of a glycine methyl ester molecule. Subsequently, the fluorescence of an N-terminal coupled fluorophore increases and can be continuously monitored (excitation wavelength 330 nm; emission wavelength 418 nm). The assay principle is shown in Figure 1.

2. Intended use

Specific and sensitive measurement of TISSUE TRANSGLUTAMINASE (TG2) activity in biological samples like tissue homogenates. Suitable for determination of TG2 concentrations of 1.25 to 10 nM (100 to 780 ng/mL or 6.4 to 110 U/L *) using a REFERENCE TG2 preparation.

* For the definition of Unit, please refer to chapter 10

3. Test sample

Sample containing human TISSUE TRANSGLUTAMINASE (hTG2). The assay is also suitable for the determination of mouse tissue transglutaminase (mTG2) and rat tissue transglutaminase (rTG2) activity.

Please note: In case of measurement of undiluted biological samples, optionally, consider the addition of dithiothreitol (DTT) to the REAGENT MIXTURE (**RM**). For details, please refer to section 7.




Figure 1: Cleavage of a covalent bond in the side chain of a specific peptidic backbone is catalyzed by tissue transglutaminase. The release of the dark quencher is followed by incorporation of glycine methyl ester. After the dark quencher is released, the fluorescence of the fluorophore increases.

The assay continuously monitors TG2 activity.

4. Reagents in the kit

- (1) A163/SR SUBSTRATE REAGENT (**SR**):
2 x 23 µL TG2 specific assay peptide (DMSO solution)
- (2) B146/BR BUFFER REAGENT (**BR**):
2 x 18 mL TRIS buffer pH 7.5 containing 10 mM calcium chloride,
100 mM sodium chloride, 5 mM glycine methyl ester (lyophilizate)
- (3) D146/DB DILUTION BUFFER (**DB**): 1 x 18 mL 10 mM TRIS buffer pH 8.1 containing 150
mM sodium chloride, 1 mM EDTA, 5 mM DTT (lyophilizate)

Optionally available at Zedira (not included in the kit):

- T022 hTG2 REFERENCE, Human tissue transglutaminase
- T040 mTG2 REFERENCE, Mouse tissue transglutaminase
- T038 rTG2 REFERENCE, Rat tissue transglutaminase
- A168 ABZ-PEPTIDE-CALIBRATOR,
reaction product formed for calibration purposes
- Z006 “Z-DON”, TG2-inhibitor

5. Reagent preparation, storage, and stability

ZEDI*X*CITE 330/418 FLUOROGENIC TG2-ASSAY KIT F014 (BUFFER REAGENT **BR** and DILUTION BUFFER **DB** lyophilizates and SUBSTRATE REAGENT **SR**) must be stored at 2-8°C. Shipment at ambient temperature is possible. The unopened reagents are stable according to the retest date printed on the box.

Table 1: Reconstitution of Kit components

Component	Preparation	Storage
BUFFER REAGENT (BR) 2 x 18 mL TRIS buffer lyophilizate (50 mM TRIS, 10 mM CaCl ₂ , 100 mM NaCl, 5 mM Gly-OMe, pH 7.5)	Add 18 mL of deionized water per vial and mix carefully	Consume within one day or store frozen at -20°C for at least 4 weeks
DILUTION BUFFER (DB) 1 x 18 mL TRIS buffer lyophilizate (10 mM TRIS, 150 mM NaCl, 1 mM EDTA, 5 nM DTT, pH 8.1)	Add 18 mL of deionized water per vial and mix carefully	Consume within one day or store frozen at -20°C for at least 4 weeks
SUBSTRATE REAGENT (SR) 2 x 23 µL TG2 specific assay peptide (DMSO solution)	Ready-to-use DMSO solution Bring to ambient temperature	Store at 2-8°C Stable for at least 12 months

Table 2: Preparation of REAGENT MIXTURE (RM)

Component	Preparation	Storage
SUBSTRATE REAGENT (SR) BUFFER REAGENT (BR)	Add 20 µL SR to 18 mL prepared BR and mix Optionally, add 90 µL of DTT (1 M in water, details: section 7)	Consume within 2 hours Protect from light Keep RM at ambient temperature until usage

Table 3: Components optionally available at Zedira (not included in the kit)

Component	Preparation	Storage
T022 hTG2 REFERENCE	Dissolve as indicated in the respective CoA Serial dilution 100-6.25 nM (7.80-0.05 µg/mL) in DB Final assay concentration will be one-tenth	Prepare enzyme diluent directly before use Do not vortex! Store diluent on ice
T040 mTG2 REFERENCE	Dissolve as indicated in the respective CoA Serial dilution 100-6.25 nM (7.80-0.05 µg/mL) in DB Final assay concentration will be one-tenth	Prepare enzyme diluent directly before use Do not vortex! Store diluent on ice
T038 rTG2 REFERENCE	Dissolve as indicated in the respective CoA Serial dilution 100-6.25 nM (7.80-0.05 µg/mL) in DB Final assay concentration will be one-tenth	Prepare enzyme diluent directly before use Do not vortex! Store diluent on ice
A168 ABZ-PEPTIDE-CALIBRATOR reaction product formed for calibration purposes	Prepare a DMSO stock solution or serial dilution according to your desired assay concentration as indicated in the respective PDS	DMSO stock solutions can be stored at -20°C for at least 6 months
Z006 TG2-INHIBITOR "Z-DON"	Prepare a DMSO stock solution or serial dilution according to your desired assay concentration as indicated in the respective PDS	DMSO stock solutions can be stored at -20°C for at least 6 months



6. Schematic assay overview

Reagent preparation

Step 1




Step 2




Sample
(not included)

Assay procedure (96 well microtiter plate)



Assay procedure (cuvette)



7. Procedure and Equipment

Set the fluorescence spectrophotometer temperature to 37°C, if applicable. Data shown in the manual are obtained at this temperature. However, the kit may also be run at ambient temperature.

The ZEDIXCITE 330/418 FLUOROGENIC TG2-ASSAY KIT (F014) can be used in fluorescence plate readers using microplates as well as in standard fluorescence spectrophotometers with cuvettes. Refer to the instructions of the manufacturer.

Add SAMPLE (**S**) and REAGENT MIXTURE (**RM**) depending on your assay format:

Microtiter plate (96 well, 300 µL): Select a microplate that is rated for fluorescence-based assays and exhibits little or no autofluorescence in the emission range of the reagent you wish to use. Black plates are typically recommended.

Preload the wells with 30 µL of your SAMPLE (**S**). Start the reaction by adding 270 µL of REAGENT MIXTURE (**RM**) to SAMPLE (**S**) and mix thoroughly.

Fluorescence cuvette (1 mL): Start the reaction by adding 100 µL of SAMPLE (**S**) to 900 µL of REAGENT MIXTURE (**RM** - prewarmed to 37°C), mix thoroughly.

Optionally: In case of measuring undiluted biological samples, consider the addition of dithiothreitol (DTT) to the REAGENT MIXTURE (**RM**). Briefly, add 90 µL of DTT (1 M in water) to 18 mL of **RM** and mix thoroughly.

Blank: Use DILUTION BUFFER (**DB**) instead of SAMPLE (**S**) to generate a blank. Measurement of samples in duplicate or triplicate is recommended.

Start the kinetic measurement 3 minutes after starting the reaction using the instrument settings shown in table 4.

Table 4: Instrument settings for fluorescence spectrophotometer.

Excitation wavelength	330 nm
Emission wavelength	418 nm
Assay time (min)	30

Determine the slope of fluorescence increase over a reaction time of 30 min. Use the linear part of slope for assessment.

Unlike absorbance, fluorescence is not an absolute measurement. The RFU scale cannot be standardized. Accordingly, absolute counts cannot be compared between readers of different manufacturers. The intensity of a fluorescent signal is usually relative to other measurements, to a reference measurement, or to the gain settings [3].

Adjust fluorescence gain of the sample with the expected highest signal output (e.g. your positive control) to avoid saturation but still cover a good assay dynamic range.

NOTE: A TG2 REFERENCE is optionally available (see section 4) for quantification. You may also consider using Zedira product A168 (“ABZ-PEPTIDE-CALIBRATOR”), the reaction product formed, for calibration purposes.

8. Number of sample measurements

The assay reagents per kit are sufficient for 130 measurements in microtiter plates and 38 measurements in cuvettes.

9. Results

Fluorescence increase is proportional (2nd polynomial fit) to the TISSUE TRANSGLUTAMINASE activity. The correlating TG2 amount can be calculated with a reference curve. Figures 2 to 4 show typical plots.




Figure 2A: **Recombinant human Tissue Transglutaminase** (T022, hTG2 REFERENCE) dependent increase in fluorescence emission (RFU: relative fluorescence units) for a serial dilution from 780 to 50 ng/mL (10 to 0.63 nM) over 25 min. The increase in fluorescence was found to be linear between 10 and 20 min. Each concentration was determined in triplicate. The control without hTG2 did not yield any increase in signal intensity (not shown). The graph represents the average curves of triplicate measurements (error bars are not shown due to readability).

Figure 2B: Plotting Δ RFU/min (taken from figure 2A, assessed between 10 and 20 min) against the indicated hTG2 concentration, a non-linear regression fit (2nd polynomial fit, $R^2 > 0.99$) was obtained. The lower limit of quantitation of the assay was found to be 100 ng/mL (1.25 nM) hTG2.




Figure 3: Reference curve for **recombinant mouse TISSUE TRANSGlutAMINASE (T040, mTG2 REFERENCE)**, diluted with DILUTION BUFFER (DB).




Figure 4: Reference curve for **recombinant rat TISSUE TRANSGlutAMINASE (T038, rTG2 REFERENCE)**, diluted with DILUTION BUFFER (DB).

10. Calibration

Since the determination of a hTG2 concentration (in ng/mL or nM) might be not applicable, the enzyme activity can be calculated to **Enzyme Units** using Zedira product A168 (“ABZ-PEPTIDE-CALIBRATOR”).

The standard definition for 1 unit (U) is the amount of enzyme that catalyzes the reaction of 1 μmol of substrate per minute. However, in most applications, the conversion of 1 μmol of substrate is not feasible and other definitions may be preferred.

In this assay the common non-standard definition of **1 U = 1 nmol/min** is used.

The conversion of enzyme concentration to units requires the extinction coefficient, determined by the measurement of the fluorescence (RFU) of a serial dilution of A168 (“ABZ-PEPTIDE-CALIBRATOR”) as shown in figure 5. The micromolar extinction coefficient ϵ $[\mu\text{M}]^{300\mu\text{L}}$, measured with 300 μL assay volume in a microtiter plate, equals to the slope of the linear regression and was determined to

$$\epsilon [\mu\text{M}]^{300 \mu\text{L}} = 743,995 \left[\frac{\text{RFU} \times \text{mL}}{\text{nmol}} \right]$$

Note: Only valid for the individual Fluorescence Microplate Reader; cannot be translated to any other reader. Determination of ϵ is required on each individual reader.




Figure 5: Determination of the extinction coefficient ϵ $[\mu\text{M}]^{300\mu\text{L}}$ by plotting the fluorescence (RFU) of a serial dilution of A168 (“ABZ-PEPTIDE-CALIBRATOR”). The extinction coefficient equals to the slope of the linear regression.

In figure 2B, the *y*-axis ($\Delta\text{RFU}/\text{min}$) represents the increase of fluorescence corresponding to the reaction product formed, plotted against the indicated hTG2 concentration.

The *x*-axis in figure 2B is converted to **Enzyme Activity** by dividing the increase of fluorescence ($\Delta\text{RFU}/\text{min}$) by the extinction coefficient resulting in figure 6.

$$\text{Enzyme Activity} = x = \frac{y}{\epsilon^{300} \mu\text{L}} \times DF \quad \left[\frac{\text{RFU} \times \text{nmol}}{\text{min} \times \text{RFU} \times \text{mL}} = \frac{\text{nmol}}{\text{min} \times \text{mL}} = \frac{\text{U}}{\text{mL}} \right]$$

$$\text{with } y = \frac{\Delta\text{RFU}}{t} \quad \text{and} \quad \text{dilution factor } DF = \frac{V}{v} = \frac{\text{total assay volume}}{\text{volume of sample}}$$

In this example, e.g. 8,177 $\Delta\text{RFU}/\text{min}$ correspond to an **Enzyme Activity** of 0.11 U/mL or 110 U/L. This also corresponds to 10 nM hTG2.

However, in case there is no reference TRANSGLUTAMINASE available, measurement of the extinction coefficient (figure 5) is sufficient to determine the **Enzyme Activity** of biological samples like tissue extracts.




Figure 6: Enzyme Activity of a hTG2 REFERENCE in U/L, determined by dividing the increase of fluorescence ($\Delta\text{RFU}/\text{min}$) by the extinction coefficient.

NOTE: All examples shown are only valid for the individual Fluorescence Microplate Reader used for measurement and cannot be translated to any other reader.

Reference measurements or calibration is required on each individual reader. Suitable references and a calibrator are optionally available at Zedira, but not included in the kit.

11. Reference Range

The ZEDI*X*CITE 330/418 FLUOROGENIC TG2-ASSAY KIT (F014) is suitable for measurements of 1.25 to 10 nM (100 to 780 ng/mL or 6.4 to 110 U/L) of TG2.

Please note: ZEDI*X*CITE 330/418 FLUOROGENIC TG2-ASSAY KIT is able to detect - but not to quantify - TG2 concentrations lower than 100 ng/mL in e.g. biological samples.

12. Kinetic parameters

The determination of TISSUE TRANSGLUTAMINASE activity is based on the ability of TG2 to cleave isopeptide bonds. The kinetic parameters obtained for TG2 specific peptidic substrate in ZEDI*X*CITE 330/418 FLUOROGENIC TG2-ASSAY are summarized in table 5. For comparison, the non-specific fluorogenic substrate A101 [2] is presented. In summary, the ZEDI*X*CITE substrate peptide performs about 7-fold better than A101.

Table 5: Kinetic parameters of TG2 substrates determined by ZEDI*X*CITE 330/418 FLUOROGENIC TG2-ASSAY with hTG2.

	K_M [μ M]	k_{cat} [s^{-1}]	k_{cat}/K_M ($M^{-1} s^{-1}$)
ZEDI <i>X</i> CITE substrate peptide A163	5.7	0.01	1,804
A101	8.6	0.002	252

13. Limitations

The ZEDI*X*CITE 330/418 FLUOROGENIC TG2-ASSAY KIT (F014) is meant for research and development only. The kit has been optimized for the measurement of purified TG2 in buffer. In case of determination of TG2 in biological samples, suitable controls are recommended (e.g. spiking with purified TG2 as positive control or blocking TG2 with Z006, a TG2 inhibitor, as negative control).

14.Precision

The coefficient of variation determined during kit development was 3.4 % for 5 nM (390 ng/mL) recombinant human TISSUE TRANSGLUTAMINASE (hTG2). Day to day variance was 1.1 %.

References

[1] Parameswaran, K.N. *et al. J. Biol. Chem.* **1997**, 272, 10311.

[2] Oertel, K. *et al. Anal. Biochem.* **2007**, 367, 152.

[3] BMG Labtech “How to optimise fluorescence gain”