

# Product Data Sheet



Product number **Various (see table below)**

Revision number **RN2.1**

**Product Name** Cereal Protein Extracts, for product numbers and names see table below.

	Barley	Rye	Wheat	Wheat Durum	Spelt	Oat	Millet (Sorghum)	Rice	Corn	Soy*
<b>Albumin + globulin</b>	G018	G021	G036 Leukosin + Edestin	G039	G030	G033	G045	G027	G024	G042
<b>Prolamin</b>	G019 Hordein	G022 Secalin	G037 Gliadin	G040	G031	G034 Avenin	G046 Kafirin	G028 Oryzin	G025 Zein	G043
<b>Glutelin</b>	G020 Hordein	G023 Secalinin	G038 Glutenin	G041	G032	G035 Avenalin	G047	G029 Oryzenin	G026 Zeanin	G044

\* corresponding protein extracts

**Background info** Cereals have a protein content of about 10%, which are classified by their solubility according to Thomas Burr Osborne (1919, The vegetable proteins.):

- Albumins + globulins: soluble in saline
- Prolamins: soluble in ethanol
- Glutelins: soluble in propanol/urea/DTE

Prolamins and glutelins like gliadin and glutenin from wheat play a crucial role in celiac disease pathophysiology.

**Description** Cereal proteins were extracted according to the method of Wieser *et al.* as published in Ruh *et al.*, 2014.

Albumins and globulins are lyophilized from 0.4 M NaCl, 67 mM Na<sub>2</sub>HPO<sub>4</sub>, 67 mM K<sub>2</sub>HPO<sub>4</sub>, pH7.6.

Prolamins are lyophilized from 60% ethanol.

Glutelins are lyophilized from 50% Propanol-1 (v/v), 2 M Urea, 10 mM DTE, 50 mM Tris HCl, pH 7.5.




Figure 1:  
Silver stained SDS-PAGE gel of **albumin + globulin** fractions (5 µg protein per lane). Protein content was determined with Bradford Assay.

\* Corresponding protein extracts

Product number **Various (see table below)**

Revision number **RN2.1**




Figure 2:  
Silver stained SDS-PAGE gel of **prolamin** fractions (5 µg protein per lane).  
Due to poor stainability 50 µg oat prolamin (G034), 20 µg rice prolamin (G028) and 20 µg corresponding soy protein extracts (G043) have been loaded.  
Protein content was determined by weighting the freeze-dried material.  
\* Corresponding protein extracts




Figure 3:  
Silver stained SDS-PAGE gel of **glutelin** fractions (5 µg protein per lane).  
Protein content was determined by comparison with weighted prolamins on silver stained SDS-PAGE (figure 2).  
\* Corresponding protein extracts

**Quantity** 5 mg

Protein content of albumins and globulins was determined using Bradford Assay against BSA as reference.

Prolamin quantities were determined by weighting the freeze-dried material (± 0.1 mg per vial).

Glutelin quantities were determined by comparison with weighted prolamins on silver stained SDS-PAGE.

**Appearance** White lyophilized solid

**Reconstitution** Add the volume of water or buffer specified in the certificate of analysis under aliquotation or the desired volume to the vial of lyophilized powder. Rotate vial gently until solid dissolves. After reconstitution, the solution should be stored frozen in working aliquots.

**Storage** Store at -20°C, desiccate

**Reference** Ruh *et al.*, J. Agric. Food Chem. 2014, 62:7604-11;  
Wieser *et al.*, Cereal Chem. 1998, 75:644-50

# Product Data Sheet



Product number **Various (see table below)**

Revision number RN2.1

---

**Related products**

- A011 - Monoclonal antibody to gliadin (clone XGY1)
- A057 - Monoclonal antibody to deamidated Gliadin
- A062 - Monoclonal antibody to deamidated and non-deamidated Gliadin
- A035 - Monoclonal antibodies to gliadin (Set No 1 comprising 12 gliadin antibodies: clone XGY1; XGY2; XGY4; XGY5; XGY8; XGY10; XGY12; XGY15; XGY16; XGY17; XGY23 and XGY24)

**Release date** 01 November 2023

**NOTE** INTENDED FOR RESEARCH USE ONLY, NOT FOR USE IN HUMAN, THERAPEUTIC OR DIAGNOSTIC APPLICATIONS.