

Microbial transglutaminase is widely used in food industry. Furthermore there are growing numbers of diagnostic and medical applications like cross-linking, site specific protein labeling or PEGylation.

Zedira's kit Z009 is a convenient and reliable product to determine the amount of transglutaminase activity.

- | | | |
|--|--------------|--------------|
| • Z009 Microbial Transglutaminase Assay Kit | 1 kit | 295 € |
|--|--------------|--------------|

Related products:

- | | | |
|--|--------|-------|
| • T001 Microbial transglutaminase, highly purified | 25 U | 255 € |
| • C100 MTG-Blocker (inhibits microbial transglutaminase activity) | 25 mg | 365 € |
| • A019 Polyclonal antibody to microbial transglutaminase | 500 µg | 345 € |



Already got our new **catalogue 2012?**

Download your electronic copy
or send us an e-mail to receive
your personal hard copy.

Get more information about transglutaminases and Zedira on Facebook  and Twitter  !

For detailed information concerning our products please refer to our web site: www.zedira.com.

Should you have any questions or queries about our products, please do not hesitate to contact us: contact@zedira.com.

We are constantly enlarging our product portfolio in order to fulfill researchers' requirements in the field of transglutaminase.

The newsletter is sent to known interested parties only. If you do not want to receive the Transglutaminase Newsletter anymore, please click here to [unsubscribe](#) your mail address.

In case that you know potential new interested colleagues, please forward this e-mail.

This newsletter is published by:

Zedira GmbH

Roesslerstr. 83

64293 Darmstadt

Germany

Phone: +49 6151 3251-00

Fax: +49 6151 3251-19

Web: www.zedira.com

E-mail: contact@zedira.com

© 2012, Zedira GmbH